## **Service Manual**

**Heraeus Multifuge X3, X3F, X3R, X3FR**

**Heraeus Megafuge 40, 40R,**

**Sorvall Legend XT, XF, XTR, XFR**

**Sorvall ST 40, ST 40R**

**Thermo Scientific SL40, SL40R**

The original version of this manual is valid. Thermo Scientific cannot guarantee the accuracy or completeness of subsequent copies

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## TABLE OF CONTENTS

Section	Title	Doc.-No.	Page
<b>1</b>	<b>OPERATING INSTRUCTIONS</b>		
<b>2</b>	<b>SERVICE</b>		
2.1	Servicing Schedule		2-1/3
2.2	Trouble Shooting		2-4/5
2.3	Error Code		2-6/10
2.4	Test Points		2-11
2.5	Imbalance Behavior		2-12
2.6	Cleaning of Instrument Parts		2-12
2.7	Electrical Safety Check		2-13
<b>3</b>	<b>FUNCTIONAL DESCRIPTION</b>		
3.1	Block Functions		3-1
3.2	Main Board Functions		3-1/3
3.3	Indication Board		3-3
3.4	Speed Detection Board, Imbalance sensor		3-4
3.5	Lid lock system		3-4
3.6	Cooling Plant		3-5
<b>4</b>	<b>DIAGRAMS</b>		
4.1	Terminal Connection		4-1
4.2	Wiring Diagram Air-Cooled		4-2
4.3	Wiring Diagram Refrigerated		4-3
4.4	Main Board, Component Plan		4-4
<b>5</b>	<b>DISASSEMBLY OF INSTRUMENT PARTS</b>		
5.1	Housing / Casing Parts		5-1/2
5.2	Electrical Components		5-3/4
5.3	Drive Components		5-5
5.4	Components of Cooling Plant		5-6
<b>6</b>	<b>SPARE PART FIGURES AND LISTS</b>		
6.1	Spare Part Figures – Air-Cooled		6-1/13
6.2	Spare Part Lists - Refrigerated		6-14/28
<b>7</b>	<b>TECHNICAL BULLETINS</b>		
<b>8</b>	<b>Preventive maintenance checklist</b>		
	<b>Calibration Certificate</b>		

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## 2.1 Servicing Schedule (yearly procedure recommended)

### 2.1.1 Maintenance Routine without Dismantling the Centrifuge

#### 2.1.1.1 Electrical Installation and Safety

- Switch OFF the centrifuge and disconnect the unit from power, check voltage supply and mains fusing (16 Amps, slow blow characteristic)
- Check condition of plug and wall socket - (let) replace defective parts
- Check cord condition and fixing / connection - replace or refit it
- Check condition of instrument socket and replace it in case of bad contacts

#### 2.1.1.2 Location and Mechanical Installation

- Check the base (ground, table, lorry with lockable wheels etc.) for resonance-free and stable conditions
- Check for a well ventilated place and sufficient distances to walls or adjacent equipment, without exposition to direct sunlight
- Check the levelling of the centrifuge drive with use of a spirit level

#### 2.1.1.3 Lid Tumbler Mechanism and Safety Device

- Connect the centrifuge to power and switch ON
- Check for correct automatic lid closing and self-acting lid opening - if in disorder, readjust lid, swivel hinges, and/or locking devices
- Check the rubber gasket for lid's and bolt's sealing and replace, if damaged
- For checking the electronically safety circuit: start the centrifuge, let it shortly run and stop it, the lid must not be unlocked by the microprocessor as long as the speed is more than 60 rpm - if safety circuit is out of function, replace main board

#### 2.1.1.4 Cleanliness of Spin Chamber and Motor Casing

- Open the lid and remove the rotor by pressing the AutoLock
- Clean the spin chamber with a dry and absorbent cloth (remove all dust and moisture - see also section for Cleaning)
- Check the cleanliness of the motor flange and take care of the annular slot around the motor shaft: penetrating fluids can damage the upper spindle bearing or electronics, remove fluids with an injector and/or absorbent paper

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**Servicing Schedule****2.1.1.5 Rotor and Accessories Condition and Sealing**

- Check the condition of rotors and accessory parts (especially all supporting or stressed partitions): the rotor and/or accessory parts must not be used any longer, if there are visible traces of mechanical damage or rust
- Check the condition of rotor and/or accessory sealing and replace them in case of malfunction

**2.1.1.6 Rotor Fixing and Drive Spindle**

- Check the perfect condition of the AutoLock system in case of malfunction order the Tool Kit # **70904693** and change the AutoLock.
- Check the condition of the drive cone: the centrifuge must not be used any longer, if the drive cone is damaged (e.g. bend or it's bearings are worn out)

**2.1.1.7 Temperature Level**

- Check air inlet slots for free ventilation, insufficient air flow will lead to an inadmissible temperature rise of compressor, motor and electronic parts

**2.1.1.8 Imbalance Behaviour**

- install available and empty rotors and check the imbalance behaviour with rotor dependant cut off and run through weights (see 2.4) and replace worn out motor rubber mounts, faulty main board or sensor board

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**Servicing Schedule**

## 2.1.2 Maintenance Routine after Dismantling the Centrifuge Casing

2.1.2.1 *Motor Supporting Elements*

- Check the supporting and damping elements of the drive motor once a year and replace them in case of increased rubber abrasion or abundance of imbalance but at least every 3 years

2.1.2.2 *Braking Circuit*

- Check the function of the braking circuit (warming up of brake resistor, even and noiseless brake effect) and replace defective parts in case of malfunction
- The break resistor is only changed by changing the complete electronic cassette

2.1.2.3 *Lead and Screwing Connections*

- Check the terminal and plug connections of all leads and on all boards and electrical components, tighten all loosen screwing connections, refit or replace defective parts
- Check the wedged connections of the board, mechanical and electrical components and re-tighten them if necessary (use screw locking lacquer for motor mounts and lid lock assembly)

2.1.2.4 *Protection Earth Core and Grounding Connections*

- Check the protection earth core for continuity and all grounding plug connectors (see also section 2.7 Electrical Safety Check)
- Check isolation resistance and accessible current (see section 2.7)

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2.2 Trouble Shooting

Error Indication	Error Cause	Possible Error Source	Corrective Procedure
Displays remain dark	no mains voltage supply	mains switch with integrated fuse or fuses on main board	switch on again, disconnect electric. components success., search for other faults
		mains fuse or circuit breaker failed	check fuse or circuit breaker and replace or switch on again
		faulty mains cord or instrument socket	check instrument cord and socket, replace defective parts
	no low voltage supply for indication board	faulty connection from CPU to indication board	check connections on CPU, indication board and connecting leads, replace defective parts
faulty indication or CPU board		replace main board completely	
display shows <b>BOOTING</b>	CPU program reset may be caused by EMI	reduced voltage supply (>10%)	remedy the failure if the voltage drops often, use a voltage stabiliser
		bad or missing ground connection	check all ground connections and the ground connection of all boards
drive makes <b>noises</b> -no good separation result	mechanics	wear out of motor rubber mount	replace motor rubber mounts (at least every three years)
	electrical	defective terminal connection, faulty lead or motor winding	winding resistances -see test points
		faulty electronic	replace main board
“PUSH OPEN”	locking drive is not supplied with voltage	faulty driving circuit	check voltage for drive motor check switch. after power ON replace the main board
	Latch motor	faulty motor winding	replace complete latch
	lid is not correctly locked	lid bolt is jamming	push lid centrally into lock and press the key again
lid is deformed		re-adjust the lid centrally	

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## 2.3 Error-Code

Error Indication	Error Cause	Possible Error Source	Corrective Procedure
„E-01“ message appears in speed display	System clock pulse not stable	hardware fault on main board	change main board
„E-02“ message	Program sequence was disturbed	internal Software problem	replace main board
„E-03“ message	Stack overflow	internal Software problem	no action in the field possible
„E-04“ message	Stack underflow	internal Software problem	no action in the field possible
„E-05“ message	Reset by HW watchdog timer	software problem	no action in the field possible
„E-08“ message	Check sum error E-PROM	E-PROM failed	replace main board
„E-09“ message	Check sum error NV RAM	NV RAM failed	replace NV RAM, or main board
„E-10“ message	NV RAM not initialised	initialising of NV-RAM is not valid	check NV-RAM and socket, insert the correct NV-RAM
„E-11“ message	NV RAM doesn't match E-PROM	check versions of NV RAM and E-PROM	replace NV RAM
„E-12“ message	Check sum error NV RAM	NV RAM failed	switch OFF/ON, replace NV RAM
„E-13“ message	Check sum error NV RAM	NV RAM failed	replace NV RAM
„E-14“ message	Chamber over temperature: display >50°C	compressor didn't work	check leads, terminals, start relay, capacitor, cooling system
		faulty cooling cycle	exp. valve, refrigerant loss
„E-15“ message	Open temp. sensor	check sensor, control circuits	replace sensor or main board
„E-16“ message	Short cut temp. sensor	check sensor, control circuits	replace sensor or main board

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Error Indication	Error Cause	Possible Error Source	Corrective Procedure
“E-44”	DC link voltage too low	line voltage <10% faulty mainboard	correct line voltage replace mainboard
“E-45”	Safety test during restart failed	during spin, switch OFF/ON, speed detection fault	change speed detection, motor plug detached
“E-61”	Latch motor not disconnected	relay not disconnect latch motor	replace main board
„E-70“	CLASS A trap software break	wrong parameter in the controller	replace main board
„E-71“	Software reset	internal controller failure	replace main board
„E-72“	Incorrect rotor recognition	software failure	no action in the field
„E-73“	Incorrect user interface identification	user interface not known by NV-RAM	replace NV-RAM
„E-74“	Motor control circuit failed	control circuit not set	replace main board
„E-75“	Faulty RAM	faulty data content in the controller	replace main board
„E-76“	Wrong information DC bus	relay not disconnect motor after lift lid	replace main board, check overtemp. switch of the resistor
„E-77“	Fault in configuration check	software faulty configured the controller	change main board
„E-78“	Set speed > max speed	signal for speed fault	replace main board

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## 2.4 Test Points

Test Points	Unit Value	Conditions
Voltage input board's resistance	120/230V 390Ω	all given values are related on 230V ( $\pm 10\%$ ) mains voltage for board's supply
motor winding resistance 20°C -insulation value	3x2,2Ω > 10MΩ	switch OFF unit, pull off motor plugs, measure inter U -> V, V -> W, W -> U resistance inter each phase and motor casing <b>Caution!</b> Motor has no PE
plugs <b>XC15/16</b> lid motor	5Ω	condition at rest, resistance at 20°C
Compressor terminal <b>X4</b> terminal box -insulation value	230V: 3,9 Ω 14,6 Ω  120V: 1,5 Ω 7,4 Ω  > 10MΩ	Switch OFF unit, compressor temp. 20°C Resistance of main coil + (starting relay 0.1Ω) Resistance of auxiliary coil (joints 13 -> 14) Resistance inter phase and motor casing
Fan cooling plant terminal <b>X2</b>	100 Ω 82 Ω	Switch OFF unit, fan resistance temperature 20°C wire black to blue wire brown to blue
Temp. measuring plug <b>X11</b>	Depends on tempera- ture	Chamber temperature 20°C (PT1000: $\approx 1078\Omega$ ) Chamber temperature -10°C (PT1000: $\approx 961\Omega$ ) Chamber temperature 40°C (PT1000: $\approx 1155\Omega$ )
speed detection  speed:	Pulses at LED Rear panel	board with active filter circuit low voltage supply light barrier is passed through (slot position) light barrier is out (turn cone a little by hand) when centrifuge is running

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## 2.5 Imbalance Behaviour

- **Install** in table available **rotor(s)** in unloaded condition

Rotor	Run through weight	Cut off weight
75003607/75003608	20 grams	35 – 45 grams
75003617	20 g	35 – 45 g
75003602	4 g	6 – 7 g
75003661	5 g	10 – 15 g

- Insert rotor dependent **cut-off weights** in all buckets or Sample holders one after the other and the centrifuge must cut off with “Imbalance” message
- Insert the **admissible imbalance weights** in the same manner, the centrifuge must run through 4 times to maximum rotor speed

## 2.6 Cleaning of Instrument Parts

### ATTENTION - WARNING!

**The electrical and electronic components must not be cleaned with moist detergents!**

For Cleaning the centrifuge housing or its accessories see Operating Instructions section 1 (maintenance and care).

- **Electronic components**  
Clean dusty components carefully with a dry and soft brush and remove loose dust with a vacuum cleaner.
- **Fans**  
Scratch off carefully with a knife or similar tool the crusted dirt from the fans' blades. Resulting grooves or marks must be removed subsequently with a fine abrasive cloth. Loose dirt is to remove with an absorbent cloth or vacuum cleaner.
- **Vent holes**  
Remove dirt from the vent grid behind the front panel by using a brush and vacuum cleaner.

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## 2.7 Electrical Safety Check

### ATTENTION!

A final electrical safety check must be performed after each maintenance and/or repair!

- **Resistance check of protective conductor**

The measuring value of the resistance between the mains plug's grounding pin and the grounding conductors of the motor, electronic chassis and the casing must not exceed 200 mΩ


Caution! The motor doesn't have PE connection. The red wire is only potential equalisation.

- **Insulation resistance Check**

Check also the insulation resistance between the poles of the mains plug and the grounding conductor; the resistance value must be more than 2 MΩ

- **Accessible current measured to IVD medical equipment regulation**

The accessible current must not exceed **0.5 mAmps** in single fault condition (interrupted protection earth wire)!

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### 3 Functional Description

#### 3.1 Block Functions

The Multifuge X3/X3R, Legend XT/XTR is a refrigerated/non-refrigerated microprocessor controlled laboratory tabletop centrifuge with noise reduced induction drive motor incorporating automatic rotor and imbalance detection and self-closing door locking systems.

The unit is equipped with following boards and components (block diagram 4-1):

- Instrument socket for mains cable at the back (units for 200V or 208V mains supply are equipped with step-up transformer)
  - Two poles mains switch with integrated thermal over-current release F1 (16A)\*
  - Main board with microprocessor part and power electronics
  - indication board with backlight liquid crystal display and a film keyboard
  - Sensor board with rotor detection sensors
  - 3 phase induction motor with integrated thermal overtemp. switch (C. O. 150°C)
  - Speed detection and imbalance board with filter attached to the motor's bottom by screws
  - self-closing door locking system (mechanical bolt keeper, motor dis-locking)
  - One break resistor mounted on the electronic cassette
  - Cooling plant with R404 compressor, expansion valve and over-pressure switch\*
  - axial-flow fans for power electronics, motor and cooling plant
- \*Multifuge X3R and Legend XTR only

#### 3.2 Main Board Functions

The **Main board** is mounted on a aluminum cassette on the rear of the centrifuge behind the rear panel. The components on main board are arranged in following function blocks

- Noise filter in accordance with EN 55011 for motor and triac circuits
- Power pack for low voltage supply of microprocessor part and power electronics
- Triac driving circuits for compressor cut-in relay, for 2 independent fan wiring terminals and for connection of a heater (special use)
- DC intermediate circuit for supply of frequency converter equipped with automatic loading and voltage (brake) controlling circuits
- Power part inclusive physically separated driving stages
- Microprocessor part with CPU and NV RAM
- Driving component for self-closing lid lock device
- Exchangeable NV-Ram containing specific data of the unit

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## Main Board Functions

### 3.2.1 Power Pack

The power pack consists of 2 transformers, 2 rectifiers (single diode) and 2 voltage regulators, which generates the following low voltage levels:

- U1 = 5V: By switched-mode regulator A8499, VDD, VCC  
supplies all processors and its electronic components on main and indication board
- U2 = 24V: Unregulated,  
is used for the driving of lid lock motor and relays K1, K3;  
the reference potential is GND, VSS
- U3 = 20V: By 2nd transformer, rectifier diodes and fix voltage regulator 7815 serves for supply of driving of the power electronics via the protection circuit by closed lid switches  
reference potential: **0V CAUTION - MAINS POTENTIAL!**

### 3.2.2 Intermediate Circuit

The DC intermediate circuit serves as an energy store between the AC power input and the transmitted motor performance. The intermediate circuit consists of a heat-sink-cooled bridge rectifier and several reservoir capacitors which are softly charged across a power resistor R2 when ever the unit is switched on. After the loading time (CPU initialisation time) the power resistor is linked by the relay contacts of K2. When switching the unit off, the DC voltage of the intermediate circuit is discharged across a path of resistors.

### 3.2.3 Brake Path

Electrical power is fed back into the intermediate circuit during motor deceleration (motor acts as generator). This braking power is transformed into heat by the connected resistor so that the intermediate circuit voltage does not rise to an excessive level. The brake resistor is switched into the intermediate circuit by a fast switching transistor (pulse-width modulation, synchronised by the double mains frequency). This transistor is voltage dependent controlled by a self-acting stage (closed loop).

### 3.2.4 Frequency Converter

The motor is 3-phased, provided with chopped direct voltage blocks (chopped frequency approx. 5kHz, amplitude is height of DC intermediate voltage). These blocks are variable in frequency and pulse-width modulation, 120° out of phase. These 3 phases are controlled during acceleration, running at set speed and deceleration (e.g. for small speed: low frequency and small pulse-width length will be affected). The FC is protected against over current, over temperature and over-voltage.

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## Main Board Functions

### 3.2.5 Microcontroller (Central Processing Unit) Part

System Menu: press any key during switch on to enter the system menu. The system menu has the submenus: language, end beep, keypad beep, lcd power save, (only HIGH END: lid pass code, program pass code), auto lid open, software and cycles counter.

Service Menu: to enter service menu press program key “1 and 2” during switch on. The submenus are: cycles (rotor cycles), error list (last 10 errors), statistic (operating hour meter, motor hour meter, cooling system hour meter, lid lock cycles).

ROUTINE: “SPEED” key to move in the menu, “BUCKET” to choice, “START” key to confirm.

HIGH END: “ENTER” to get to the next submenu, numeric to choice. For “YES” it’s “1” for “NO” it’s “0”. To confirm press “START” key.

The actually indicated version No.'s are displayed in the time field and will be upgraded if program or data changes are necessary.

The control program (ROM) is integrated in the controller component. The most important operating parameters (e.g. the maximum data for all rotors, the last operator settings, the basic values for speed and temperature controlling, data of acceleration and deceleration curves as well as data sets of program places are stored in the non volatile (NV)-RAM.

### 3.2.6 Temperature Measuring

The temperature inside of the rotor chamber is measured by a PT1000 sensor. The sensor's head is mounted in the bowl. The signal is led across XU to the PCB where a pulsed current by Q1 through 0,1R - 100R - 0,1R and PT1000 given to the CPU. The switched input voltage is internally converted into data for further use.

### 3.3 Indication Board

The indication board is mounted behind the operating panel. The connection to the main board is done by a RJ45 cable. Data transfer as well the low voltage supply of 5V is verified. The operating keys are managed by an own connection.

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### 3.4 Speed Detection Board, Imbalance sensor

On a small board, which is attached in a cap on the motor, the circuits for speed detection, imbalance and rotor recognition are arranged. A flashing LED at low speed shows the function of the circuit. A measurement of the signal is not possible.

In the rotor's bottom 4 magnets are forced, being alternately arranged in polarity and positioned to each other to different distances on a constant circular ring (segment steps of 20°). A Hall sensor is exactly fixed below to the sensor board. The magnets are switching a flip-flop stage by which different pulse lengths and pulse intermissions are generated during one rotor rotation. By this signal the CPU is able to detect up to 65 different rotors and the correct direction of rotation, too.

Additional in the cap is a acceleration detector fixed. Dependent of the rotor imbalance this sensor generates an analogous voltage signal which is led via a following band-pass filter to main board, where it is analysed by CPU. When the CPU has identified the installed rotor type, the programmed speed dependent imbalance limitation curve is loaded. The analogue imbalance signal is filtered, converted into digital form and continuously compared with this specific limitation curve. If the actual imbalance signal will rise at any speed beyond the limit, the error signal for imbalance load "imbalance" will be displayed and the unit slows down.

### 3.5 Lid lock system

The latch is directly controlled by CPU via a PWM motor driver. If all necessary signals are existing which allow an opening the motor is fired with a pulsed voltage of 24 VDC. For additional safety the relay K11/12 will be activated if the speed is below 120 rpm controlled by triggered mono-flops.

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
### 3.6 Cooling plant (only Multifuge X3R, Legend XTR)


The Multifuge X3R / XTR is equipped with a high-grade cooling plant, to meet its requirements of driving power and compact design. The CFC free refrigerant R404A evaporates at low temperatures ( $\leq -20^{\circ}\text{C}$ ) and absorbs through this the frictional heat from chamber produced by the fast spinning rotor. 5 fans serve for even air flow through condenser so that the refrigerant's condensing temperature is kept on a low level (lateral front fans are only switched on, if  $n > 0$  rpm).


A specially designed coaxial heat exchanger serves for essential undercooling of the condensed refrigerant to ensure an optimum evaporation and a perfect operation of the expansion valve. To maintain low evaporation and resultant condensation temperatures on one hand and to support the starting procedure of the compressor on other hand, a thermostatic expansion valve with 10bara MOP for 230V Centrifuges and 5,0bara MOP for 120V Centrifuges is used.


The minimum achievable sample temperatures are depending on speed, rotor type and a little on ambient temperature. After a compensating time, given by the pre-cooling or pre-heating program as optimal result, the rotor has reached the selected set temperature within a tolerance of  $\pm 2\text{K}$ .

To reduce the power consumption of the unit the compressor will not be activated during accelerating of the motor. When reaching the selected speed the compressor will be activated if a low temperature had been selected.

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Due to cooled and non-cooled version are described not all parts can be available.

## 5.1 Dismantling the Housing

### 5.1.1 Front Panel

- Open the lid, switch off the unit and **pull out the mains plug**
- remove the two screws on top of the front-panel
- lift the front-panel and take it off
- Disconnect wiring from display board

### 5.1.2 Lid

- Open the lid, switch off the unit and **pull out the mains plug**
- Remove rear-panel
- Remove gas pressure spring from top holder – see 5.1.5
- Remove the 2 screws for lid's attaching on the casing and the 4 from below
- Take the screw driver out of joint blocks to separate lid from gas stay and remove the lid from casing
- When reassembling the lid take care that both bolt keepers engage the lock entries concentric and stand out of the screw-in hollow exactly for 46mm

### 5.1.3 Casing

- Remove the front panel - see 5.1.1
- Remove the lid - see 5.1.2
- Remove the 3 remaining screws on both sides and one from both lid lock.
- remove the casing with lid and bolt holes sealing in vertical position from armoured chamber and frame
- Reassemble the casing analogously in reverse order and take care of the correct lid sealing seat

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#### 5.1.4. Automatic Lid Locking System

- Dismantle the front panel - see 5.1.1
- Loosen the 4 screws from above
- Disconnect plug from the motor
- Reinstall the latch and gear motor in reverse order

#### 5.1.5 Gas lid stay

- Open the lid, switch off the unit, **pull out the mains plug**
- Remove rear-panel
- Knock carefully the joint bolts for upper attaching of gas lid stay out of joint block **with its three notches ahead** by use of a hammer and pin punches (but press the bolt in with non notches ahead by use of pliers)
- Push an adequate screw driver into bolt's place to keep the gas pressure spring in position to avoid lid's slamming
- Remove screws at the bottom trough upper notch
- Take out the gas pressure spring to the top
- Knock carefully the joint bolt for lower attaching of gas lid stay out of joint block **with its three notches ahead** by use of hammer and pin punches (but press the bolt in with non notches ahead by use of pliers)
- Replace worn out gas lid stay and reinstall the new one analogously in reverse order and remount the unit completely.

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## 5.2 Replacement of Electrical Components

### 5.2.1 NV-RAM on the Main Board

- Notice old displayed identification and version numbers of NV-RAM when powering on
- Remove main board - see 5.2.2
- **ATTENTION - CMOS components! Discharge your body before handling!** Pull NV-RAM carefully out of socket using chip removal tool
- Reinsert the new IC correctly
- Put the electronic cassette back in position, connect motor with the electronic assy switch power ON, check the displayed identification and version numbers of the new ICs and refit rear-panel and the motor flange (pay attention of position from temperature sensor only X3R/XTR)

### 5.2.2 Main Board Replacement

- Dismantling the motor flange and deposit it in the chamber (don't remove)
- Disconnect the motor cable of the electronic cassette by catching the motor cable between the motor and chamber (if there is not enough space remove the motor – see 5.3.1)
- turn the centrifuge and remove the rear-panel
- Remove all wires from the terminals by pulling of
- Loose the 4 screws which holds the aluminum cassette with the electronic and pull it out
- **Touch a grounded receptacle to discharge your body before touching the sensitive CMOS components!** Take the NV-RAM out of box (or re-use the old one when in **perfect** condition) and insert it into socket of the new main board (if necessary, remove the placed but non-programmed NV-RAM before)
- Remount new electronic cassette
- **Do not mix up disconnected cables during re-connection!**
- Reassemble the device in reverse order and perform a test run, **making sure the drive turns in the right direction** (see imprinted arrow on rim of rotor chamber)!
- Watch correct position of the ferrite core rings

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### 5.2.3 Key and Display Board

- Remove the front panel - see 5.1.1
- Remove the attaching screws and take out the indication board
- Unplug the connecting line to main board and transfer to new board
- Reassemble the indication board analogously in reverse order, check controller identification and refit the front panel.

### 5.2.4 Sensor Board and Speed Detection Board

- The sensor board is placed with the Speed Detection Board under the cover on the motor and fixed with 3 screws.
- Loose the main board to disconnect the cables – see 5.2.2
- Remove auto lock system – **see instruction 50119249**
- Reassemble the sensor board analogously in reverse order

### 5.2.5 Speed Detection Board

- See – 5.2.4

### 5.2.6 Fans

- **Condenser fan, right side:**
- Remove front panel - see 5.1.1
- Remove the casing – see 5.1.3
- Remove the screws of fan to condenser
- Take fan apart

### 5.2.7 Brake Resistor – not available as single component

- Remove the electronic cassette – see 5.2.2 and install a new one
- break resistor can only be changed by installing a new electronic cassette
- Check the cleanliness of air louvers and clean them if necessary, pay attention to good heat transmission
- Reassemble the device analogously in reverse order. Perform some short test runs with brake curve 9

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### 5.3 Replacement of Drive Components

#### 5.3.1 Drive Motor

- Unscrew motor take it off
- Disconnect the wires for motor from electronic cassette
- Remove auto-lock system from motor – **see instruction 50119249**
- Remove Motor cover with recognition cover
- Re-install recognition cover to new motor

#### 5.3.2 Motor rubber Mounts

- Remove motor – see 5.3.1
- Loosen with a 13mm fixed spanner the screws under the AV mount, by positioning the fixed spanner between the ventilation slot of the base plate. Don't forget the 4<sup>th</sup> AV mount at the side
- Inspect the AV mounts and if necessary replace them
- All 4 AV mounts must be exchanged completely at once!

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**5.4 Replacement of Parts of the Cooling Plant**

First of all is to figure out that high ambient temperatures ( $\vartheta_U > 35^\circ\text{C}$ ), insufficient heat transfer (e.g. polluted condenser, de-heater, air gaps) or electrical errors (like defective fans, sensor or main board, starting relay or capacitor) can be excluded from the cooling problem.

Pressure / temperature ranges are valid for continuous cooling at max. speed,  $\vartheta_U = 25^\circ\text{C}$ !

Procedure for replacing defective components of the cooling cycle:

- Dismantle the casing -see 5.1.3 and mount pierce valves onto the service lines
- When gauging too little static pressure perform a leakage test ( $\text{N}_2$  max. 25bars)
- Discharge the refrigerant with suction device into recycling bottle
- Cooling plant is bonded by Exp.valve with evaporator and compressor suction side with evaporator
- Now loose the 6 screws which fixed the cooling plant with the base plate
- Remove the hole cooling plant
- In case of compressor fault, let the oil run out into a bowl to recycle it
- Unsolder defective
- Let a little amount of dry nitrogen gas stream through the tubes when soldering to avoid scaling on the tube's inner surface
- Reassemble the cooling plant
- Do the evacuation and refill the system with the exact amount of liquid refrigerant on the pressure side via a new external dryer (check amount with a balance or filling cylinder, re-solder the service lines tight after finishing work)
- Reassemble the device and perform a cooling power test run

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Charging information 230 V version	Charging information 120 V version
Refrigerant: R 404a Amount: 520 grams	Refrigerant: R 404a Amount: 520 grams
Operating pressure with rotor 75003607 and round buckets 75003608	Operating pressure with rotor 75003607 and round buckets 75003608
Speed: 4700 rpm Ambient: 23 °C Time: 3 hours	Speed: 4700 rpm Ambient. 23°C Time: 3 hours
Suction: about 3,6 bars (-13°C) Discharge: about 18 bars (38°C)	Suction: about 2,5 bars (-13°C) Discharge: about 9,5 bars (30°C)

*The world leader in serving science*

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Heraeus Megafuge 40 / Sorvall ST 40

Heraeus Multifuge X3 / Sorvall Legend XT

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


Gas spring P/N 20310534

[Pos. 4000](#)

*The world leader in serving science*

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Heraeus: Multifuge X3R, Megafuge 40R  
Sorvall: Legend XTR, ST 40R

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- [Outer Casing – Aussengehäuse](#) [Parts List](#)
- [Lid Cover – Deckel](#) [Parts List](#)
- [Electronics – Elektronik](#) [Parts List](#)
- [Service Kit](#) [Parts List](#)
- [Cooling System – Kälteanlage](#) [Parts List](#)

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# Calibration Certification Procedure

## Multifuge 3X and versions

Revised - 04/12/2008

### 1. General

Successfully completing the Preventive Maintenance Checklist and Calibration Certification Procedure as outlined will ensure that the instrument is properly maintained and calibrated to Thermo Scientific specifications.

### 2. Equipment Required

- Digital Multimeter
- Temperature Meter w/probe
- Stopwatch
- Photo Tachometer
- Safety Tester
- Service Manual
- Centrifuge Operating Manual

**Note:** Check calibration due date on all test equipment before starting certification.

### 3. Documentation Required

- Preventive Maintenance Procedure and Checklist (Service Manual)
- Calibration Certification Procedure and Form
- Certification Label
- Certificate of calibration for each piece of test equipment used to perform the certification
- Training Certificate

#### Preventive Maintenance Check

- Perform Preventive Maintenance checks as outlined in the Service Manual to ensure the instrument is in good working order without performing any calibrations.
- Complete the Preventive Maintenance Checklist. (calibration data will be filled in after Certification)

#### Calibration Certification

- **SPEED**

Install a rotor, set the speed control to 3000 rpm (or wanted) and start the centrifuge. When the speed has stabilized, measure the actual speed with an external tachometer (photo tachometer on multimeter where speed = frequency[no conversion]). Record measured speed and displayed speed in the speed section of the certification form. Set the speed control to max. rpm and after stabilization, record measured speed and displayed speed in the speed section of the certification form. Determine if the speeds recorded meet the specifications on the form. If the specifications are met, proceed to the next section. A recalibration of the instrument is not possible.

# **Calibration Certification Procedure**

## **Multifuge 3X and versions**

- **TEMPERATURE**

Place water into the rotor or bucket. Select a temperature of +4°C, maximum allowed speed of the rotor and time of minimum 3 hours. After that record the displayed temperature, measured temperature inside rotor or bucket after rotor has come to rest, in the temperature section of the certification form. Repeat the test at 15°C. Record the values in the temperature section of the certification form. Determine if the temperature values recorded meet the specifications on the form. If the specifications are met, proceed to the next section. A recalibration of the instrument is not possible. Initial the temperature section. To pre-cool the rotor select pre-cooling program of the centrifuge.

- **TIME**

Using a stopwatch, check the timer by setting a 10 minute run. Start the centrifuge and measure from the time the start button is pushed until the timer switches to stop. Record the measured time in the time section of the form. Determine if the time value recorded meets the specifications in the form. If the specifications are met, proceed to the next section. A recalibration of the instrument is not possible. Initial the time section.

- **SAFETY TEST**

Using safety tester and connect the centrifuge to it. Measure Insulation resistance.. Select to Earth Conductor resistance and touch with probe all mechanical part of the centrifuge. Resistance must be <0.2 ohms at all points otherwise check for good grounding. Measure accessible current to procedure in the service manual.

### **Documentation**

After filling in the certification information complete the following documentation:

- Log test equipment data on calibration certification form.
- Review P.M.checklist and calibration certification form to ensure that all data and required information has been recorded properly.(including your signature and date)
- Fill out a certification label and place it on the right side of the centrifuge where visible or at a location selected by the customer.
- Ask the customer to sign the calibration certification form and place a copy in the customer's Certification Record Book for that centrifuge.
- Place copies of the P.M.checklist and Field Service Report in the same Record Book.
- Provide copies of your training certificate and test equipment calibration certificates if requested by the customer.

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# Thermo Scientific Service Information

**Kommentar:**

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Datum \_\_\_\_\_

# Thermo Scientific Service Information

## Heraeus Centrifuges PREVENTIVE MAINTENANCE CHECKLIST

<b>Multifuge X3 &amp;versions</b> <b>Type</b> <b>S/N</b>
--

### Rotor Inspection

- o Check for corrosion and wear
- o Autolock
- o Inspect rotor magnets (rotor I.D.)  
Serial numbers:  
75003607.....  
75003668 .....
- 75003670.....
- 75003603.....
- 75003606.....
- 7500.....

### Preliminary Checks

- o Clean inst. & condenser coils
- o Rotor chamber (defrost)
- o Inst. environment  
(level & clearance)

### Pre-run Checks

- o Ground continuity
- o Line voltage \_\_\_\_\_V AC
- o Display on power-up
- o Door interlock and latch
- o Drive stop & Door unlatch (<50 rpm)
- o Fan operation (condenser & drive)

### Mechanical

- o Motor bearings
- o Microprocessor Initialization (if required)
- o Refrigerant leaks (if required)
- o Compressor current draw

### Performance Checks

- o Acceleration time \_\_\_\_\_seconds
- o instrument top speed \_\_\_\_\_rpm
- o Deceleration time \_\_\_\_\_seconds
- o Time (10 minute set) \_\_\_\_\_actual min:sec
- o Temperature (4° set) \_\_\_\_\_actual (°C)

### Imbalance

- o Centrifuge stable and level
- o Centrifuge shuts off w/proper weight

### Safety test

- o Accessible leakage current <500µA
- o Insulation resistance >10 Megaohm
- o Earth Conductor Resistance <0.2 Ohm

# Thermo Scientific Service Information

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SPEED / DREHZAHL						
Standard Test			Rotor Used:			
Set Speed (rpm)	Tolerance (whichever is greater)	Displayed Speed (rpm)	Measured Speed (rpm)	Meets Spec. (y/n)	Notes	Passed (Initials)
4000	10 rpm					
15000	10 rpm					
Custom Test			Rotor Used:			
	10 rpm					
	10 rpm					

TEMPERATURE / TEMPERATUR						
Standard Test			Rotor Used:			
Calibration Temp. (°C)	Tolerance (K)	Displayed Temp. (°C)	Measured Temp. (°C)	Meets Spec. (y/n)	Notes	Passed (Initials)
+4	+/- 2					
+15	+/- 3					
Custom Test			Rotor Used:			
	+/- 2					
	+/- 2					
	+/- 2					

TIME / ZEIT						
Standard Test			Rotor Used:			
Timer Set (min.)	Tolerance (sec.)	Measured Time (m/s)	Meets Spec. (y/n)	Notes	Passed (Initials)	
10	+/- 30					
Custom Test			Rotor Used:			
	+/- 30					
	+/- 30					
	+/- 30					

TEST EQUIPMENT / PRÜFMITTEL			
Type	Model	Serial Number	Date of Next Calibration
Stopwatch			
Temperature Meter			
Photo Tachometer			
Safety Tester			

SAFETY TEST / SICHERHEITSTEST		
Insulations resistance	Earth Conductor Resistance	Accessible Current
>10 Megaohm	<0.2 Ohm	<500µA

**Revised 04/12/08**  
 Certified By: \_\_\_\_\_ Reviewed By: \_\_\_\_\_ Expiration Date \_\_\_\_\_

Signature \_\_\_\_\_ Date \_\_\_\_\_ Customer Signature \_\_\_\_\_ Date \_\_\_\_\_  
 Use ONLY Black permanent ink when filling out. Leave NO blank spaces.